



WESTERN INTERNATIONAL BAND CLINIC

AUDITION PACKET - SET A

GETTING STARTED:

1. Students need to record all three pieces from SET A:
 - Classical Solo - varies by instrument
 - Chromatic Scale
 - Danny Boy
2. Use any recording equipment software or app to record the auditions:
 - Mac - Quicktime, GarageBand or Audacity
 - PC - Quicktime, Sound Recorder or Audacity
 - Smartphone/Tablet - Search the app store for “free MP3 recorder”
3. Export the auditions in MP3 or M4a format.
4. Create an account at the www.bandworld.org/wibc website or login
5. Upload all three files from SET A audition music.

NOTE TO PERCUSSIONISTS: You must record TWO of the three percussion auditions (snare, mallets and/or timpani). Mallet players do not have to record “Danny Boy” in addition to the normal mallet audition music found at the back of this packet.

ALL INSTRUMENTS: It is okay to use All-Northwest or your All-State Music instead of WIBC audition music.

NOTE TO DIRECTORS: There is another set of audition materials, in addition to audition packet SET A. If you would like to rehearse all of the WIBC audition pieces at the same time, with your full band, download the SET B audition materials from our website and use the scores and parts provided. See the info below on SET B audition packet for more info.

(Download Audition Packet - SET B at www.bandworld.org/wibc - click on “Honor Band Audition Information” left-hand column.) Audition Packet - **SET B** was created with the full band in mind during rehearsals. You can handout the individual parts and play them all together during class time. To allow for optimal instrument ranges, the various audition pieces are in two keys at the same time. The harmonies are perfect fourths or perfect fifths to also foster additional tuning possibilities during full band rehearsal. Our hope is that you will utilize these audition pieces during warm-up time throughout the school year. Your band students will become familiar with the scale, minor etude and “Danny Boy” which will prepare them to audition for WIBC. We hope you find these new sets of WIBC audition materials useful. Please feel free to share the PDFs or our website with other band directors. You also have permission to use these materials for your own auditions.

All auditions must be uploaded to www.bandworld.org/wibc

QUESTIONS? Contact Scott McKee • (541) 778-4880
scottmckee@bandworld.org • www.bandworld.org/wibc

WESTERN INTERNATIONAL BAND CLINIC AUDITION

(SET A)

Allegro in C Minor

Bass Clarinet

G. F. Handel

Allegro ♩ = 84 - 96

mf

8

f

12

f

16

19

21

2

25

p

mf

pp

26

mf

29

30

f

32

3

36

f

p

38

cresc.

40

poco a poco

41

4

47

f

poco rit.

WESTERN INTERNATIONAL BAND CLINIC AUDITION

(SET A)

Bass Clarinet

Chromatic Scale

$\text{♩} = 96$

5

The Chromatic Scale is written in treble clef with a 4/4 time signature. It consists of two lines of music. The first line starts with a quarter note G4, followed by eighth notes G#4, A4, A#4, B4, B#4, C5, C#5, D5, D#5, E5, E#5, F5, F#5, G5. The second line starts with a quarter note G5, followed by eighth notes G#5, A5, A#5, B5, B#5, C6, C#6, D6, D#6, E6, E#6, F6, F#6, G6. The piece ends with a whole note G6.

Danny Boy

Old Irish Air

Slowly and Freely

5

9

13

a tempo *rit.* *Slower* *rit.*

The Danny Boy piece is written in treble clef with a 4/4 time signature and a key signature of two sharps (D major). It consists of four lines of music. The first line starts with a quarter note D4, followed by eighth notes D4, E4, F4, G4, A4, B4, C5, B4, A4, G4, F4, E4, D4. The second line continues with eighth notes C4, B3, A3, G3, F3, E3, D3, C3, B2, A2, G2, F2, E2, D2. The third line starts with a quarter note D3, followed by eighth notes D3, E3, F3, G3, A3, B3, C4, B3, A3, G3, F3, E3, D3. The fourth line continues with eighth notes C3, B2, A2, G2, F2, E2, D2, C2, B1, A1, G1, F1, E1, D1. The piece ends with a whole note D1. Performance markings include *Slowly and Freely* at the beginning, *a tempo* at measure 13, *rit.* at measure 14, *Slower* at measure 15, and *rit.* at measure 16.